

ERK1 / ERK2 Monoclonal Antibody

Catalog No.	A10613	Category	Monoclonal Antibodies
Applications	WB, IHC	Observed MW	39kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	36kDa/41kDa/38kDa/40kDa/43kDa

Immunogen Information

Immunogen	A synthetic peptide of human ERK1 / ERK2
Gene ID	5594/5595
Swiss prot	P28482/P27361
Synonyms	

Product information

Source	Mouse
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:1000 -
1:2000

IHC 1:50 -
1:200